

Instructions for Use

ROTRON N



Tem Innovations GmbH

tem®

Intended Use:

ROTRON N is a quality control material for monitoring accuracy and precision of tests carried out on the ROTEM® delta Thromboelastometry System. For *in vitro* diagnostic use only.

Reagents:

Product Name:

Reference Number:

Package size:

5 x 1 vial ROTRON N Lyo

5 x 1 vial ROTRON N Dil (1300 µl)

After reconstitution, each vial is sufficient for 4 assays.

Constituents:

ROTRON N control plasma (lyophilized) was made from human plasma collected using sodium citrate as an anticoagulant (3.2%). The plasma was adjusted to yield coagulation values in the ROTEM® approximating the normal range of whole blood (Level 1 control). Stabilizers and buffers were added prior to the lyophilization.

Preparation of the ready-to-use control:

Let the ROTRON N Dil and Lyo vials reach room temperature. Dissolve the contents of the ROTRON N Lyo vial by pouring the contents of the ROTRON N Dil vial into the lyophilisate. A small drop of the diluent remains in the vial. Do not attempt to transfer the diluent with a pipette! Do not use any other diluent than the one supplied. Close the vial with the rubber cap and the screw cap and swirl gently. Take care that the powder is completely dissolved. Then let stand for 10 - 15 minutes in the closed container, to allow the plasma to reconstitute. Before use, bring to 37 °C on the temperature controlled work area for 5 minutes and carefully mix again by swirling gently.

Storage and Stability:

Store at +2 to +8°C (35.6 - 46.4 °F). The unopened ROTRON N reagent is stable until the expiry date indicated on the label.

Stability after Initial Use:

Reconstituted ROTRON N is stable for 8 hours at 2 - 8 °C (35.6 - 46.4 °F). Avoid contamination and always close the vials again (rubber and screw cap) after each use. Freezing and thawing of ROTRON N is not recommended.

Additional Material:

ROTEM® device; Cup & Pin (measurement cells; REF 200011); pipette tips (REF 400040 / REF 400041); in-tem® reagent (REF 503-02-US) for intrinsic activated thromboelastometry, ex-tem® reagent (REF 503-05-US) for extrinsic activated thromboelastometry, and star-tem® reagent (REF 503-10-US) for recalcification. The reagents are used in accordance with the instructions in the respective package insert.

Method:

Analytical Principle:

The use of controls is a valuable tool to ensure the quality of the coagulation tests. The package insert of ROTRON N shows batch-specific target ranges. The material can therefore not only be used to monitor precision but also the accuracy of the system (device, reagents, user) (2, 3).

Limitation of the Procedure:

All controls are subjected to the limitation of the test system. Variables such as temperature, reagent stability, instrument properties and individual techniques may affect the final result. Always strictly follow the manufacturer's guidelines for the device and reagents.

Target ranges:

Each batch of ROTRON N is supplied with a batch-specific quality certificate with a table of target ranges for the individual tests. When using ROTRON N and the Tem Innovations system reagents, the results for the specified tests should be within these ranges. If other reagents are used, reagent-specific target ranges must be generated by the user.

Measurement Calculation:

The ROTEM® device calculates numerous parameters. These parameters and their mathematical background can be found in the ROTEM® operating manual. The quality certificate for ROTRON N only gives target ranges for the important primary parameters. Most of the other parameters are derived parameters which are derived by conversion from the primary parameters.

Warnings:

Precautions:

Each donor unit used in the preparation of ROTRON N has been tested for antibodies against HIV Type 1 and 2, Hepatitis C-Virus antibodies, Treponema pallidum antibodies as well as Hepatitis B surface-antigen and Hepatitis C genome by PCR. The tests used are all CE certified tests according to list A of the European Directive for IVDs (98/79/EC) and are under supervision by the responsible European governmental authority. The plasmas were found to be negative on the tested parameters. However, since no test can completely rule out the presence of blood borne diseases, these control plasmas have to be handled as potentially infectious material (4).

For In Vitro Diagnostic Use only

Recommended External Control Procedure:

- ⇒ Conducting a quality control serves as verification of the whole system; device, pipette, reagent, user. It is recommended to run a control test (QCIN for INTEM and QCEN for EXTEM) prior to assays of patient samples. Control tests should also be run:
 1. After preventive maintenance is performed on the instrument.
 2. Whenever a new shipment or lot number of reagent is received and prior to use on the ROTEM® system.
 3. Whenever a value from the ROTEM® Service Menu is out of range (Phase shift, Variance, Center, Amplitude, or Temperature) and it cannot be corrected (see ROTEM® delta operating manual chapter 5.5.1).
- ⇒ If a result is outside the target range, the test should be repeated in accordance with the test-specific instructions on the same channel and on a further second channel.
- ⇒ If both measurements are within the target range, it is likely that there has been procedural error during the first measurement.
- ⇒ If both measurements are outside the target range, it is possible that a reagent has deteriorated or there is a systematic user error. If the control is also outside the target range with fresh reagents, a defective control is likely. If the result is repeated with a freshly reconstituted control, it is likely that there is a device failure, in which case, you should contact your technical service provider.
- ⇒ If the first result on the original channel is confirmed and the result on the second channel is OK, it can be assumed that there is a channel-specific problem or a calibration problem. Do not use this channel for any further measurements. In this event, proceed as recommended in the ROTEM® user manual or contact your technical service provider.

Recommended External Control Interval:

Manufacturer recommendation: minimum every week (e.g. QCIN on channel 1 - 2 and QCEN on channel 3 - 4 with ROTRON N the first week and QCINP and QCENP with ROTRON P the following week).

CLIA recommendation: minimum daily and/or within every 8 hours that a patient test is resulted. Additionally testing should be performed if required. External quality control should be performed at 2 levels (e.g. QCIN on channel 1, QCEN on channel 2 with ROTRON N; QCINP on channel 3, QCENP on channel 4 with ROTRON P).

CLIA equivalent QC testing, option 2 (QSA.02.04.01): Internal quality control (see ROTEM® delta operating manual chapter 5.5.1) should be checked at a minimum daily and external quality control should be run at a minimum weekly using 2 levels of controls (e.g. QCIN on channel 1, QCEN on channel 2 with ROTRON N; QCINP on channel 3, QCENP on channel 4 with ROTRON P).

Performance data:

Precision:

| INTEM / QCIN | CT CV[%] | α-angle CV[%] | A20 CV[%] |
|-------------------------------|-------------|------------------|--------------|
| Within-run ¹ | 1.4 | 0.6 | 8.1 |
| Between Operator ² | 2.7 | 0.6 | 5.4 |

| EXTEM / QCEN | CT CV[%] | α-angle CV[%] | A20 CV[%] |
|-------------------------------|-------------|------------------|--------------|
| Within-run ¹ | 3.1 | 0.6 | 11 |
| Between Operator ² | 7.9 | 0.3 | 5.5 |

- ¹ .5 runs on each of the 4 channels of one instrument using ROTRON N
- ² 5 operators run ROTRON N in duplicates

Bibliography:

- (1) NCCLS Document H21-A2. Collection, transport, and processing of blood specimens for coagulation testing and performance of coagulation assays, 3rd ed. Approved Guideline 1998
- (2) Blutgerinnungsstudien mit der Thrombelastographie, einem neuen Untersuchungsverfahren. Hartert, H.; Klijn, Wochenschrift 1948; 26: 577-583
- (3) Thromboelastographic Coagulation Monitoring during Cardiovascular Surgery with the ROTEM Coagulation Analyzer, Calatzis, A. et al.; Management of Bleeding in Cardiovascular Surgery edited by Roque Pilarre; Hanley & Belfus, Inc. Philadelphia, PA, 2000
- (4) Biosafety in Microbiological and Biomedical Laboratories, U.S. Department of Health and Human Services, Washington 1993 (HHS publication No. (CD) 93-6395)

Manufacturer:

Tem Innovations GmbH
Martin-Kollar-Str. 13-15
D-81829 Munich / Germany
Email: info@tem-innovations.de

Distributor USA:

TEM Systems, Inc.
3210 E. NC Highway 54 Suite 150
RTP, NC 27709 USA
Email: info@roteminc.com

2015-10
V 0012
EN USA

| ??? | Symbols on labels | Symboles sur les étiquettes | Simbolos en las etiquetas |
|----------------|----------------------------------------------------|------------------------------------------------------|---------------------------------------------------------------|
| LOT | lot | Lot | Lote |
| REF | reference or order number | Référence ou numéro de commande | Referencia o número de pedido |
| | store between x and y degree celsius | A conserver entre x et y °C | Conservar entre "x" y "y" grados Celcius |
| CE | CE marking according to IVD guideline 98/79 EC | Marquage CE conforme aux directives IVD 98/79 EC | Marca CE según directiva IVD 98/79 CE |
| | expiry date | Date d'expiration | Fecha de caducidad |
| CERT | certificate (standard curve and evaluation table) | Certificat (courbe standard et tableau d'évaluation) | Certificado (curva estándar y tabla de evaluación) |
| | instructions for use | Mode d'emploi | Instrucciones de uso |
| IVD | in-vitro diagnostic use | Diagnostic in vitro | Para uso en diagnóstico in vitro |
| | single use only | A usage unique | Úsese solo una vez |
| SN | serial number | Numéro de série | Número de serie |
| | manufacturing date | Date de fabrication | Fecha de fabricación |
| | name and address of the manufacturer | Nom et adresse du fabricant | Nombre y dirección del fabricante |
| 10 x 7 | number of tests (n) in the kit (e.g. 10 x 7 tests) | Quantité de tests (n) dans le kit (ex. 10 x 7 tests) | Cantidad de tests (n) contenidos en el kit (ej. 10 x 7 tests) |
| | biohazard | Risque biologique | Riesgo biológico |
| RCNS | reconstitute with... | A reconstituer avec... | Reconstituir con... |
| Reagent | name of the reagent (e.g. Aqua = distilled water) | Nom du réactif (ex. Aqua = eau distillée) | Nombre del reactivo (ej. Aqua = agua destilada) |
| 200 µl | volume of the diluent | Volume du diluent | Volumen del diluyente |